



Defence Research and  
Development Canada    Recherche et développement  
pour la défense Canada



# **Inactivation of Microorganisms by Gamma Irradiation**

*Bacillus atropeus spores and Erwinia herbicola*

R.E. Hilsen, B. Kournikakis, B. Ford  
DRDC Suffield

**DISTRIBUTION STATEMENT A**  
Approved for Public Release  
Distribution Unlimited

Technical Memorandum  
DRDC Suffield TM 2005-236  
December 2005

Canada

# **Inactivation of Microorganisms by Gamma Irradiation**

*Bacillus atrophaeus* spores and *Erwinia herbicola*

R.E. Hilsen, B. Kournikakis, B. Ford  
Defence R&D Canada – Suffield

**Defence R&D Canada – Suffield**

Technical Memorandum


DRDC Suffield TM-2005-236

December 2005

AQ F06-10-6964


**20060713096**

Author



Rayanne E. Hilsen


Approved by



Dr. L. P. Nagata

AI Head, Chemical & Biological Defence Section

Approved for release by



Dr. P. A. D'Agostino

Chair, DRDC Suffield Document Review Panel

## Abstract

---

The purpose of this study was to establish kill curves detailing the sterilization efficiency of  $^{60}\text{Co}$  gamma irradiation on non-pathogenic vegetative and spore preparations. This was done by testing the viability of the microorganisms over a range of gamma irradiation exposures. The procedures used to carry out the work will be used in future sterilization studies of pathogenic microorganisms for later use as antigen reagents. The initial use of non-pathogenic organisms provided a low-risk method for familiarization with the protocol and resolution of any potential problems. Results showed that 4.0 kGy irradiation over a 12 minute period was required for complete inactivation of *Erwinia herbicola*. The inactivation of *Bacillus atropheus* (formerly *Bacillus subtilis* var *globigii* (BG)) dry powder required 25kGy over 75 minutes and the inactivation of *Bacillus atropheus* in liquid suspension required 35 kGy over 105 minutes. These results showed that vegetative cells and spore preparations can be successfully sterilized by gamma irradiation in a short amount of time. These results will be used as starting points to establish kill curves for Risk Group 3 organisms such as *Bacillus anthracis* and *Yersinia pestis*.

## Résumé

---

L'objectif de cette étude consistait à établir des courbes de neutralisation détaillant l'efficacité de la stérilisation par irradiation gamma  $^{60}\text{Co}$  sur des préparations végétatives non pathogéniques et de spores. Ceci a été réalisé en testant la viabilité des microorganismes dans toute une gamme d'irradiations par expositions. Les procédures utilisées pour effectuer ces travaux serviront plus tard à l'étude de la stérilisation des microorganismes pathogéniques qui seront utilisés plus tard comme réactifs d'antigènes. L'utilisation initiale d'organismes non pathogéniques est une méthode à faible risque permettant la familiarisation avec le protocole et la résolution des problèmes potentiels. Les résultats indiquent qu'il faut une irradiation de 4.0 kGy durant une période de 12 minutes pour inactiver complètement *Erwinia herbicola*. L'inactivation du *Bacillus atropheus* (autrefois appelé *Bacillus subtilis* var *globigii* (BG) en poudre sèche exigeait 25 kGy durant une période de 75 minutes et l'inactivation du *Bacillus atropheus* en suspension liquide exigeait 35 kGy durant une période de 105 minutes. Ces résultats indiquaient qu'il est possible de stériliser les préparations de cellules végétatives et de spores par des irradiations gamma en peu de temps. Ces résultats seront utilisés comme points de départ pour établir les courbes de neutralisation d'organismes du Groupe du risque 3 tels que *Bacillus anthracis* et *Yersinia pestis*.

This page intentionally left blank.

## Executive summary

---

### Introduction

Detection assays in the laboratory often employ inactivated preparations of vegetative bacteria or spores as sources of antigens. Chemical (e.g. formaldehyde) and thermal (e.g. autoclave) inactivation can be used in the preparation of antigens but the antigen structure itself can be adversely affected by the inactivation process. The use of gamma irradiation provides better preservation of the antigenic structure while still providing a reliable method of inactivation. An MDS Nordion Gammacell 220 Excel with a  $^{60}\text{Co}$  source was used for this study.

The purpose of this study was to establish kill curves detailing the sterilization efficiency of  $^{60}\text{Co}$  gamma irradiation on non-pathogenic vegetative and spore preparations. This was done by testing the viability of the microorganisms over a range of gamma irradiation exposures. The procedures used to carry out the work will be used in future in the sterilization of pathogenic microorganisms for use as antigen reagents. The initial use of non-pathogenic organisms provided a low-risk method for familiarization with the protocol and its' refinement if problems were identified.

### Results

Results showed that 4.0 kGy irradiation over a 12 minute period was required for complete inactivation of *Erwinia herbicola*. The inactivation of *Bacillus atropheus* (formerly *Bacillus subtilis* var *globigii* (BG)) dry powder required 25 kGy over 75 minutes and the inactivation of *Bacillus atropheus* in liquid suspension required 35 kGy over 105 minutes. These results showed that vegetative cells and spore preparations can be successfully sterilized by gamma irradiation in a short amount of time.

### Significance

The gamma inactivation data in this study will be used to define the initial exposure range for further study of gamma irradiation bacterial Risk Group 3 agents (e.g. *Bacillus anthracis* and *Yersinia pestis* held by DRDC Suffield). Organisms inactivated by gamma irradiation could then be safely used as killed antigen preparations. This provides an important research tool as this type of material is not available from known commercial sources.

### Future Considerations

Since genetic repair and metabolic recovery of microorganisms after irradiation has been documented [3], sterility counts should be repeated again days later to check for any surviving organisms. It would also be beneficial for understanding the efficiency of radiation sterilization to determine the relationship between the dose of gamma irradiation delivered and the amount of organism being irradiated.

R.E. Hilsen, B. Kournikakis, B. Ford. 2005. Inactivation of Microorganisms by Gamma Irradiation – Preliminary BSL-2 Spore and Vegetative Organisms.  
DRDC Suffield TM 2005-236. Defence R&D Canada – Suffield.

# Sommaire

---

## Introduction

Les biotests de détection en laboratoire emploient souvent des préparations inactivées de bactéries végétatives ou de spores comme sources d'antigènes. L'inactivation chimique (ex : formaldéhyde) et thermique (ex : autoclave) peut être utilisée dans la préparation des antigènes mais la structure d'antigène elle-même peut être affectée de façon négative par le processus d'inactivation. L'utilisation d'irradiation gamma préserve mieux la structure antigénique tout en étant une méthode fiable d'inactivation. On a utilisé *MDS Nordion Gammacell 220 Excel* avec une source de  $^{60}\text{Co}$  pour cette étude.

L'objectif de cette étude consistait à établir des courbes de neutralisation détaillant l'efficacité de la stérilisation par irradiation gamma  $^{60}\text{Co}$  sur des préparations végétatives non pathogéniques et de spores. Ceci a été réalisé en testant la viabilité des microorganismes dans toute une gamme d'irradiations par expositions. Les procédures utilisées pour effectuer ces travaux seront utilisées dans le futur pour l'étude de la stérilisation des microorganismes pathogéniques devant être utilisés plus tard comme réactifs d'antigènes. L'utilisation initiale d'organismes non pathogéniques est une méthode à faible risque permettant la familiarisation avec le protocole et ses raffinements quand des problèmes sont identifiés.

## Résultats

Les résultats indiquent qu'il faut une irradiation de 4.0 kGy durant une période de 12 minutes pour inactiver complètement *Erwinia herbicola*. L'inactivation du *Bacillus atropheus* (autrefois appelé *Bacillus subtilis var globigii* (BG) en poudre sèche exigeait 25 kGy durant une période de 75 minutes et l'inactivation du *Bacillus atropheus* en suspension liquide exigeait 35 kGy durant une période de 105 minutes. Ces résultats indiquaient qu'il est possible de stériliser les préparations de cellules végétatives et de spores par des irradiations gamma en peu de temps.

## La portée des résultats

Les données d'inactivation gamma de cette étude seront utilisées pour définir la gamme d'expositions initiales pour les études ultérieures sur l'irradiation gamma des agents bactériens du Groupe de risque 3 (ex : *Bacillus anthracis* et *Yersinia pestis* que possède RDDC Suffield). Des organismes inactivés par irradiation gamma pourraient ensuite être utilisés de manière sécuritaire comme préparations d'antigènes morts. Ceci est un outil de recherche très important puisque ce type de matériau n'est pas disponible chez les sources commerciales.

### **Éléments à prendre en considérations pour le futur**

La réparation génétique et le rétablissement métabolique des microorganismes après irradiation ayant été documentés [3], les comptages de stérilité devraient être répétés de nouveau quelques jours plus tard pour vérifier si des organismes ont survécu. Ceci serait aussi bénéfique pour comprendre l'efficacité de la stérilisation par radiation et déterminer le rapport entre la dose d'irradiation gamma donnée et la quantité d'organismes irradiés.

R.E. Hilsen, B. Kournikakis, B. Ford. 2005. Inactivation of Microorganisms by Gamma Irradiation – Preliminary BSL-2 Spore and Vegetative Organisms.  
DRDC Suffield TM 2005-236. R & D pour la défense Canada – Suffield.



This page intentionally left blank.

## Table of contents

---

Abstract.....	i
Executive summary .....	iii
Sommaire.....	iv
Table of contents .....	vii
List of figures .....	viii
List of tables .....	viii
1. Introduction .....	1
2. Materials and Methods .....	2
2.1 Vegetative Cell Preparation.....	2
2.2 Spore Preparations.....	2
2.3 Gamma Irradiation.....	2
2.4 Viability Assays .....	3
3. Results .....	4
3.1 Inactivation of Vegetative Cells .....	4
3.2 Inactivation of Dry Spores.....	5
3.3 Inactivation of Spores In Suspension .....	6
3.4 Sample Tube Suitability .....	6
4. Discussion .....	7
5. References .....	8
ANNEX A – Gammacell 220 Decay Graph.....	9
ANNEX B – Gammacell 220 Decay Chart .....	10
List of symbols/abbreviations/acronyms/initialisms .....	11

## List of figures

---

Figure 1: Gamma Irradiation Kill Curve for <i>Erwinia herbicola</i> .....	4
Figure 2: Gamma Irradiation Kill Curve for <i>Bacillus atrophaeus</i> powder.....	5
Figure 3: Gamma Irradiation Kill Curve for <i>Bacillus atrophaeus</i> in suspension.....	6
Figure A1: Gammacell 220 Decay Graph. ....	9

## List of tables

---

Table 1: Gammacell 220 Decay Chart. ....	10
--	----

# 1. Introduction

---

Sterilization of biological products by radiation is not a new technology, but it is undergoing a revival. Irradiation of biological products by gamma rays has several advantages, including independence of packaging, excellent reproducibility, and the ability to preserve biological function such as antigenicity, protein activity [1], or even food value [2]. Gamma irradiation kills microorganisms by direct modification (e.g. crosslinking) of proteins and nucleic acids, and by free-radical effects in free water. The mechanisms of radiation effects on biomaterials (or packaging materials) could potentially compromise the structural integrity of the packaging causing leakage or breakage, thus the balance of achieving sterilization and preserving the residual structure or function of the irradiated materials must be obtained. Like other sterilization modes, chemical or physical, irradiation of various types of preparation requires validation and standardization. Ideally, the functional properties of the preparation can be preserved, while achieving sterilization to the point where no infectious or contaminating live organisms remain.

## 2. Materials and Methods

Three different experimental preparations covering both vegetative and spore forming organisms were exposed to gamma irradiation in triplicate:

### 2.1 Vegetative Cell Preparation

Overnight cultures of *Erwinia herbicola*, strain number 8L5 from the Brooks Horticultural Research Station (Brooks, AB) (2), with a stock concentration of  $8.08 \times 10^9$  cfu/mL were grown in Trypticase Soy Broth (TSB), lot number 0038002 from Becton Dickinson (Cockeysville, MD). A stock culture of 125mL of fresh TSB was inoculated with 250  $\mu$ L of *Erwinia herbicola* grown previously in TSB. *Erwinia* stock was grown overnight at 37°C with shaking at 125 rpm. Stock culture was then aliquotted into twelve 15 mL tubes with 10 mL of stock into each tube.

### 2.2 Spore Preparations

Dry Spore Preparation:

Dry *Bacillus atrophaeus* (formerly *Bacillus subtilis* var *globigii* (BG)) stock powder was weighed out into ten 15 mL tubes. A total of 0.01 g was added to each tube. The spore stock (donated by the US Department of Defense, Dugway Proving Ground, Utah), is held in dry powder form with a stock concentration of  $\sim 1 \times 10^{11}$  cfu/gm.

Spore Preparation in Liquid Media :

0.01 g of the same dry BG stock was added to 110mL sterile double-distilled water in a 250 mL flask. The preparation was mixed vigorously and aliquotted into ten 15 mL tubes with a total of 10 mL of this stock aliquotted into each tube.

### 2.3 Gamma Irradiation

Bacterial samples were gamma irradiated in 15 mL polypropylene centrifuge tubes (Fisher Scientific, Ottawa, ON). Each tube was vortexed (with the exception of the dry BG stock powder), prior to being placed into a new unlabelled metal paint can ( $\sim 10.7$ cm w x 12.4 cm h). The paint can was then loaded into the sample chamber at the top of the Gammacell unit. The irradiation timer was set for the appropriate exposure time. Gamma irradiation was carried out with the Gammacell 220 Excel  $^{60}\text{Co}$  (manufactured by MDS Nordion) as per the "Standard Operating Procedures Gammacell 220 EXCEL".

The spore preparation samples received radiation doses from 0 to 50 kGy in increments of 4.6 kGy (15 minute intervals), whereas the vegetative cell preparation received radiation doses from 1 to 4 kGy in 0.30 kGy increments (one minute intervals). The gamma source dose rate during the course of this study was 18.21 kGy per hour or 0.303 kGy/min.  $^{60}\text{Co}$  has a half life of 5.24 years so the time required to achieve a desired level of irradiation will increase with time. A chart is provided (Appendix A) for future dose calculations

## **2.4 Viability Assays**

Stock samples and gamma irradiated samples were 10-fold serially diluted in sterile double-distilled water. Total viable counts for the appropriate dilutions of either stock or gamma irradiated samples were carried out by spread plating 250  $\mu\text{L}$  of each stock or sample onto Nutrient Agar plates (Nutrient Agar powder, lot number 213000 from Becton Dickinson, Sparks, Maryland), using the Autoplate 4000 from Spiral Biotech. All stocks and samples were plated in duplicate and incubated overnight at 37°C. Plates were counted using the CASBA 4 with Colony Image Analysis (Spiral, Biotech) and the data exported to Microsoft Excel 2003 and Sigma Plot 8.0 for graphical analysis.

### 3. Results

#### 3.1 Inactivation of Vegetative Cells

A total of twelve 10 mL samples containing *Erwinia herbicola* in TBS were gamma irradiated at one minute increments. The initial stock concentration of  $3.79 \times 10^8$  cfu/mL, showed a nine log reduction at complete inactivation with 4.0 kGy. It took 12 minutes to deliver this dose to the last sample, number 12.




Figure 1: Gamma Irradiation Kill Curve for *Erwinia herbicola*. A total of 12 *Erwinia herbicola* culture samples, each in 10 mL TSB, were gamma irradiated at one minute increments. Complete inactivation is seen with 4.0 kGy at 12 minutes. This corresponds to sample number 12. Results are the average of three experiments  $\pm 1$  standard deviation.

### 3.2 Inactivation of Dry Spores

A total of ten samples containing 0.01g of *Bacillus atropheus* dry powdered spores were gamma irradiated at 15 minute increments. The initial stock concentration of  $1.06 \times 10^7$  cfu/mL showed a seven log reduction at complete inactivation with 25 kGy. It took 75 minutes to deliver this dose to sample number five.




Figure 2: Gamma Irradiation Kill Curve for *Bacillus atropheus* powder. Ten 0.01 g *Bacillus atropheus* samples were gamma irradiated at 15 minute increments. Complete inactivation is seen with 25 kGy at 75 minutes. This corresponded to sample number five. Results are the average of three experiments  $\pm 1$  standard deviation.



### 3.3 Inactivation of Spores In Suspension

A total of ten 10 mL samples containing 0.01g of *Bacillus atrophaeus* in water were gamma irradiated at 15 minute increments. The initial stock concentration of  $8.25 \times 10^6$  cfu/mL showed a seven log reduction at complete inactivation with 35 kGy. It took 105 minutes to deliver this dose to sample number seven.




Figure 3: Gamma Irradiation Kill Curve for *Bacillus atrophaeus* in suspension. Ten 10 mL samples of 0.01 g *Bacillus atrophaeus* in water were gamma irradiated at 15 minute increments. Complete inactivation is seen with 35 kGy at 105 minutes. This corresponded to sample number seven. Results are the average of three experiments  $\pm 1$  standard deviation.

### 3.4 Sample Tube Suitability

Although some of the 15 mL polypropylene tubes used for this study were exposed to gamma irradiation for over an hour, the only result of this gamma exposure was some slight yellowing. None of the sample tubes showed any visible damage due to the exposure to gamma irradiation. The appearance of the tubes exposed for the maximum amount of time indicated that the polypropylene tubes can be safely used for gamma irradiation of both wet and dry samples for at least the 150 minutes maximum dose used during this study.

## 4. Discussion

---

Detection assays in the laboratory often employ inactivated preparations of vegetative bacteria or spores as sources of antigens. Chemical (e.g. formaldehyde) and thermal (e.g. autoclave) inactivation can be used in the preparation of antigens but the antigen structure itself can be adversely affected by the inactivation process. The use of gamma irradiation provides better preservation of the antigenic structure while still providing a reliable method of inactivation. An MDS Nordion Gammacell 220 Excel with a  $^{60}\text{Co}$  source was used for this study.

A study carried out at DSTL Porton Down [4] used gamma irradiation to inactivate *Bacillus anthracis* (Sterne and Ames strains) vegetative cells and spores. This study provides us with a preliminary range of radiation exposures. DSTL Porton Down determined that vegetative cells required 2.85 kGy and spores 41.5 kGy to achieve a reliable and reproducible 100% inactivation of the organisms. The irradiation source used at Porton Down delivered a dose of 0.9 kGy per hour so that an exposure time of 50 hours was required to kill the spores. The DRDC Suffield Gammacell source is more powerful, providing a 19.4 kGy per hour exposure. Our data suggested that 100% inactivation of spores could be achieved in 105 minutes.

Inactivation of either vegetative or spore forming organisms by gamma irradiation is shown to be a convenient and reproducible technique. Organism preparations, either in suspension or in the dry state, are successfully inactivated by gamma irradiation. Vegetative cells are very sensitive to gamma radiation whereas spore forming bacteria require a much higher dose of gamma irradiation. However, the magnitude of this inactivation depends on the type of organism, the irradiation dose and physical state of organism.

The gamma inactivation data in this study will be used to define the initial exposure range for further study of gamma irradiation of bacterial Risk Group 3 agents (e.g. *Bacillus anthracis* and *Yersinia pestis*) held by DRDC Suffield. The gamma inactivated organisms could then be safely used as killed antigen preparations. This provides an important research tool since this type of material is not available from known commercial sources.

## 5. References

---

1. Tran, H., Marlowe, K., McKenney, K., Petrosian, G., Griko, Y., Burgess, W., Drohan, W., Imboden, M., Kempf, C., Boschettie, N., Mann, D. (2004). Functional integrity of intravenous immunoglobulin following irradiation with a viricidal dose of gamma radiation. *Biologicals* 32: 94-104.
2. Rocelle, M., Clavero, J., Monk, D., Beuchat, R., Doyle, P., Brackett, R.E. (1994) Inactivation of *Escherichia coli* O157:H7, *Salmonellae*, and *Campylobacter jejuni* in Raw Ground Beef by Gamma Irradiation. *Appl. Environ. Micro.* 60: 2069-2075.
3. Clavero M.R., Monk J.D., Beuchat, L.R., Doyle, M.P., Brackett, R.E. Inactivation of *Escherichia coli* O157:H7, *salmonellae*, and *Campylobacter jejuni* in raw ground beef by gamma irradiation. *Appl Environ Microbiol.* 1994 Jun;60(6):2069-75.
4. Manchee, R.J. and Watson, S. (1991). Inactivation of Vegetative Bacteria and Spores of *Bacillus anthracis* by Gamma Irradiation. (CBDE Technical Note, No 1090).

## ANNEX A – Gammacell 220 Decay Graph




Figure A1: Gammacell 220 Decay Graph. The above graph shows the decay curve for the Gammacell 220. Date at time 0 is 12 March 2004. (Received from Dr. Barry Ford).

## ANNEX B – Gammacell 220 Decay Chart

Years from Calibration	Date	CDR kGy/hr	Dose in 10 minutes kGy	Kill Curve Correction Factor	Kill Time Correction Factor
0.00	12-Mar-04	19.45	3.24	1.00	1.00
0.50	12-Sep-05	18.21	3.04	0.94	1.07
1.00	12-Mar-05	17.06	2.84	0.88	1.14
1.50	12-Sep-06	15.97	2.66	0.82	1.22
2.00	12-Mar-06	14.95	2.49	0.77	1.30
2.50	12-Sep-07	14.00	2.33	0.72	1.39
3.00	12-Mar-07	13.11	2.19	0.67	1.48
3.50	12-Sep-08	12.28	2.05	0.63	1.58
4.00	12-Mar-08	11.50	1.92	0.59	1.69
4.50	12-Sep-09	10.76	1.79	0.55	1.81
5.00	12-Mar-09	10.08	1.68	0.52	1.93
5.50	12-Sep-10	9.44	1.57	0.49	2.06
6.00	12-Mar-10	8.84	1.47	0.45	2.20
6.50	12-Sep-11	8.27	1.38	0.43	2.35
7.00	12-Mar-11	7.75	1.29	0.40	2.51
7.50	12-Sep-12	7.26	1.21	0.37	2.68
8.00	12-Mar-12	6.79	1.13	0.35	2.86
8.50	12-Sep-13	6.36	1.06	0.33	3.06
9.00	12-Mar-13	5.96	0.99	0.31	3.27
9.50	12-Sep-14	5.58	0.93	0.29	3.49
10.00	12-Mar-14	5.22	0.87	0.27	3.72
10.50	12-Sep-15	4.89	0.82	0.25	3.98
11.00	12-Mar-15	4.58	0.76	0.24	4.25
11.50	12-Sep-16	4.29	0.71	0.22	4.54
12.00	12-Mar-16	4.01	0.67	0.21	4.84
12.50	12-Sep-17	3.76	0.63	0.19	5.17

Table 1: Gammacell 220 Decay Chart. The above chart indicates gamma irradiation dose rate change over time. Due to  $^{60}\text{Co}$  having a half life of 5.24 years, the time required to achieve the initial central dose rate (kGy/hr) at 0.00 years from calibration will increase, as shown by the "Kill Time Correction Factor". (Received from Dr. Barry Ford).

## List of symbols/abbreviations/acronyms/initialisms

---

BG	<i>Bacillus atropheus</i> (formerly <i>Bacillus subtilis var globigii</i> )
BSL	Biological Safety Level
$^{60}\text{Co}$	cobalt 60
cfu/ml	colony forming units per millilitre
°C	Degrees Celsius
g	grams
h	height
kGy	kiloGrays
kGy/min	kiloGrays per minute
ul	microlitres
ml	millilitre
min	minute
rpm	revolutions per minute
TSB	Trypticase Soy Broth
w	width

UNCLASSIFIED  
**SECURITY CLASSIFICATION OF FORM**  
(highest classification of Title, Abstract, Keywords)

<b>DOCUMENT CONTROL DATA</b> (Security classification of title, body of abstract and indexing annotation must be entered when the overall document is classified)		
<b>1. ORIGINATOR</b> (the name and address of the organization preparing the document. Organizations for who the document was prepared, e.g. Establishment sponsoring a contractor's report, or tasking agency, are entered in Section 8.)  Defence R&D Canada – Suffield PO Box 4000, Station Main Medicine Hat, AB T1A 8K6	<b>2. SECURITY CLASSIFICATION</b> (overall security classification of the document, including special warning terms if applicable)  <b>UNCLASSIFIED</b>	
<b>3. TITLE</b> (the complete document title as indicated on the title page. Its classification should be indicated by the appropriate abbreviation (S, C or U) in parentheses after the title).  Inactivation of Microorganisms by Gamma Irradiation: Bacillus atropheus spores and Erwinia herbicola (U)		
<b>4. AUTHORS</b> (Last name, first name, middle initial. If military, show rank, e.g. Doe, Maj. John E.)  Hilsen, R., Kournikakis, B., Ford, B.		
<b>5. DATE OF PUBLICATION</b> (month and year of publication of document)  December 2005	<b>6a. NO. OF PAGES</b> (total containing information, include Annexes, Appendices, etc) <b>21</b>	<b>6b. NO. OF REFS</b> (total cited in document)  <b>4</b>
<b>7. DESCRIPTIVE NOTES</b> (the category of the document, e.g. technical report, technical note or memorandum. If appropriate, enter the type of report, e.g. interim, progress, summary, annual or final. Give the inclusive dates when a specific reporting period is covered.)  Technical Memorandum		
<b>8. SPONSORING ACTIVITY</b> (the name of the department project office or laboratory sponsoring the research and development. Include the address.)		
<b>9a. PROJECT OR GRANT NO.</b> (If appropriate, the applicable research and development project or grant number under which the document was written. Please specify whether project or grant.)	<b>9b. CONTRACT NO.</b> (If appropriate, the applicable number under which the document was written.)	
<b>10a. ORIGINATOR'S DOCUMENT NUMBER</b> (the official document number by which the document is identified by the originating activity. This number must be unique to this document.)  DRDC Suffield TM 2005-236	<b>10b. OTHER DOCUMENT NOS.</b> (Any other numbers which may be assigned this document either by the originator or by the sponsor.)	
<b>11. DOCUMENT AVAILABILITY</b> (any limitations on further dissemination of the document, other than those imposed by security classification)  <input checked="" type="checkbox"/> Unlimited distribution <input type="checkbox"/> Distribution limited to defence departments and defence contractors; further distribution only as approved <input type="checkbox"/> Distribution limited to defence departments and Canadian defence contractors; further distribution only as approved <input type="checkbox"/> Distribution limited to government departments and agencies; further distribution only as approved <input type="checkbox"/> Distribution limited to defence departments; further distribution only as approved <input type="checkbox"/> Other (please specify):		
<b>12. DOCUMENT ANNOUNCEMENT</b> (any limitation to the bibliographic announcement of this document. This will normally corresponded to the Document Availability (11). However, where further distribution (beyond the audience specified in 11) is possible, a wider announcement audience may be selected).  Unlimited		

UNCLASSIFIED  
**SECURITY CLASSIFICATION OF FORM**

UNCLASSIFIED  
SECURITY CLASSIFICATION OF FORM

13. ABSTRACT (a brief and factual summary of the document. It may also appear elsewhere in the body of the document itself. It is highly desirable that the abstract of classified documents be unclassified. Each paragraph of the abstract shall begin with an indication of the security classification of the information in the paragraph (unless the document itself is unclassified) represented as (S), (C) or (U). It is not necessary to include here abstracts in both official languages unless the text is bilingual).

The purpose of this study was to establish kill curves detailing the sterilization efficiency of  $^{60}\text{Co}$  gamma irradiation on non-pathogenic vegetative and spore preparations. This was done by testing the viability of the microorganisms over a range of gamma irradiation exposures. The procedures used to carry out the work will be used in future sterilization studies of pathogenic microorganisms for later use as antigen reagents. The initial use of non-pathogenic organisms provided a low-risk method for familiarization with the protocol and resolution of any potential problems. Results showed that 4.0 kGy irradiation over a 12 minute period was required for complete inactivation of *Erwinia herbicola*. The inactivation of *Bacillus atropheus* (formerly *Bacillus subtilis* var *globigii* (BG)) dry powder required 25kGy over 75 minutes and the inactivation of *Bacillus atropheus* in liquid suspension required 35 kGy over 105 minutes. These results showed that vegetative cells and spore preparations can be successfully sterilized by gamma irradiation in a short amount of time. These results will be used as starting points to establish kill curves for Risk Group 3 organisms such as *Bacillus anthracis* and *Yersinia pestis*.

14. KEYWORDS, DESCRIPTORS or IDENTIFIERS (technically meaningful terms or short phrases that characterize a document and could be helpful in cataloguing the document. They should be selected so that no security classification is required. Identifies, such as equipment model designation, trade name, military project code name, geographic location may also be included. If possible keywords should be selected from a published thesaurus, e.g. Thesaurus of Engineering and Scientific Terms (TEST) and that thesaurus-identified. If it is not possible to select indexing terms which are Unclassified, the classification of each should be indicated as with the title.)

Cobalt-60, gamma irradiation, kill curve, *Bacillus atropheus*, *Erwinia herbicola*